

Antibodies Product Datasheet

Product Name:	Nkx-2.4 Antibody
Catalog Number:	C20583
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Nkx-2.4 . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	NKX2-4
Protein Name:	Homeobox protein Nkx-2.4
Other Name:	NKX2-4; NKX2D; Homeobox protein Nkx-2.4; Homeobox protein NK-2 homolog D
Human Swiss Prot No.:	Q9H2Z4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NKX2.4 Polyclonal Antibody detects endogenous levels of NKX2.4 protein.
Storage and Stability:	-20°C/1 year