

## Antibodies Product Datasheet

|                               |                                                                                                                                                              |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HNF-1 $\beta$ Antibody                                                                                                                                       |
| <b>Catalog Number:</b>        | C20593                                                                                                                                                       |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                            |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                    |
| <b>Dilutions:</b>             | Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.                                                                          |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                             |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                            |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                        |
| <b>Immunogen:</b>             | Synthesized peptide derived from the N-terminal region of human HNF-1 $\beta$ .                                                                              |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                      |
| <b>Gene Name:</b>             | HNF1B                                                                                                                                                        |
| <b>Protein Name:</b>          | Hepatocyte nuclear factor 1-beta                                                                                                                             |
| <b>Other Name:</b>            | HNF1B; TCF2; Hepatocyte nuclear factor 1-beta; HNF-1-beta; HNF-1B; Homeoprotein LFB3; Transcription factor 2; TCF-2; Variant hepatic nuclear factor 1; vHNF1 |
| <b>Human Swiss Prot No.:</b>  | P35680                                                                                                                                                       |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                      |
| <b>Specificity:</b>           | HNF-1 $\beta$ Polyclonal Antibody detects endogenous levels of HNF-1 $\beta$ protein.                                                                        |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                 |