

Antibodies Product Datasheet

Product Name:	CRSP34 Antibody
Catalog Number:	C20596
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CRSP34 . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	MED27
Protein Name:	Mediator of RNA polymerase II transcription subunit 27
Other Name:	MED27; CRSP34; CRSP8; Mediator of RNA polymerase II transcription subunit 27; Cofactor required for Sp1 transcriptional activation subunit 8; CRSP complex subunit 8; Mediator complex subunit 27; P37 TRAP/SMCC/PC2 subunit; Transcriptional co
Human Swiss Prot No.:	Q6P2C8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CRSP34 Polyclonal Antibody detects endogenous levels of CRSP34 protein.
Storage and Stability:	-20°C/1 year