

Antibodies Product Datasheet

Product Name:	CREB3L2 Antibody
Catalog Number:	C20601
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CREB3L2. AA range:269-318
Concentration:	1 mg/ml
Gene Name:	CREB3L2
Protein Name:	Cyclic AMP-responsive element-binding protein 3-like protein 2
Other Name:	CREB3L2; BBF2H7; Cyclic AMP-responsive element-binding protein 3-like protein 2; cAMP-responsive element-binding protein 3-like protein 2; BBF2 human homolog on chromosome 7
Human Swiss Prot No.:	Q70SY1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CREB3L2 Polyclonal Antibody detects endogenous levels of CREB3L2 protein.
Storage and Stability:	-20°C/1 year