

Antibodies Product Datasheet

Product Name:	ALY Antibody
Catalog Number:	C20613
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human THOC4. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	ALYREF
Protein Name:	THO complex subunit 4
Other Name:	ALYREF; ALY; BEF; THOC4; THO complex subunit 4; Tho4; Ally of AML-1 and LEF-1; Aly/REF export factor; Transcriptional coactivator Aly/REF; bZIP-enhancing factor BEF
Human Swiss Prot No.:	Q86V81
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALY Polyclonal Antibody detects endogenous levels of ALY protein.
Storage and Stability:	-20°C/1 year