

Antibodies Product Datasheet

Product Name:	Liprin β 1 Antibody
Catalog Number:	C20615
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LIPB1. AA range:300-349
Concentration:	1 mg/ml
Gene Name:	PPFIBP1
Protein Name:	Liprin-beta-1
Other Name:	PPFIBP1; KIAA1230; Liprin-beta-1; Protein tyrosine phosphatase receptor type f polypeptide-interacting protein-binding protein 1; PTPRF-interacting protein-binding protein 1; hSGT2
Human Swiss Prot No.:	Q86W92
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Liprin β 1 Polyclonal Antibody detects endogenous levels of Liprin β 1 protein.
Storage and Stability:	-20°C/1 year