

Antibodies Product Datasheet

Product Name:	SNAT2 Antibody
Catalog Number:	C20663
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC38A2. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	SLC38A2
Protein Name:	Sodium-coupled neutral amino acid transporter 2
Other Name:	SLC38A2; ATA2; KIAA1382; SAT2; SNAT2; Sodium-coupled neutral amino acid transporter 2; Amino acid transporter A2; Protein 40-9-1; Solute carrier family 38 member 2; System A amino acid transporter 2; System A transporter 1; System N amino a
Human Swiss Prot No.:	Q96QD8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SNAT2 Polyclonal Antibody detects endogenous levels of SNAT2 protein.
Storage and Stability:	-20°C/1 year