

Antibodies Product Datasheet

Product Name:	AMY-1 Antibody
Catalog Number:	C20670
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from AMY-1 . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	MYCBP
Protein Name:	C-Myc-binding protein
Other Name:	MYCBP; AMY1; C-Myc-binding protein; Associate of Myc 1; AMY-1
Human Swiss Prot No.:	Q99417
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AMY-1 Polyclonal Antibody detects endogenous levels of AMY-1 protein.
Storage and Stability:	-20°C/1 year