

Antibodies Product Datasheet

Product Name:	FEV Antibody
Catalog Number:	C20672
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the N-terminal region of human FEV.
Concentration:	1 mg/ml
Gene Name:	FEV
Protein Name:	Protein FEV
Other Name:	FEV; PET1; Protein FEV; Fifth Ewing variant protein; PC12 ETS domain-containing transcription factor 1; PC12 ETS factor 1; Pet-1
Human Swiss Prot No.:	Q99581
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FEV Polyclonal Antibody detects endogenous levels of FEV protein.
Storage and Stability:	-20°C/1 year