

Antibodies Product Datasheet

Product Name:	MNT Antibody
Catalog Number:	C20673
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MNT. AA range:315-364
Concentration:	1 mg/ml
Gene Name:	MNT
Protein Name:	Max-binding protein MNT
Other Name:	MNT; BHLHD3; ROX; Max-binding protein MNT; Class D basic helix-loop-helix protein 3; bHLHd3; Myc antagonist MNT; Protein ROX
Human Swiss Prot No.:	Q99583
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MNT Polyclonal Antibody detects endogenous levels of MNT protein.
Storage and Stability:	-20°C/1 year