

Antibodies Product Datasheet

Product Name:	Ndfip1 Antibody
Catalog Number:	C20684
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Ndfip1 . at AA range: 150-230
Concentration:	1 mg/ml
Gene Name:	NDFIP1
Protein Name:	NEDD4 family-interacting protein 1
Other Name:	NDFIP1; N4WBP5; PSEC0192; PSEC0223; NEDD4 family-interacting protein 1; Breast cancer-associated protein SGA-1M; NEDD4 WW domain-binding protein 5; Putative MAPK-activating protein PM13; Putative NF-kappa-B-activating protein 164; Putative Q9BT67
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ndfip1 Polyclonal Antibody detects endogenous levels of Ndfip1 protein.
Storage and Stability:	-20°C/1 year