

Antibodies Product Datasheet

Product Name:	Mindin Antibody
Catalog Number:	C20686
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SPON2. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	SPON2
Protein Name:	Spondin-2
Other Name:	SPON2; DIL1; Spondin-2; Differentially expressed in cancerous and non-cancerous lung cells 1; DIL-1; Mindin
Human Swiss Prot No.:	Q9BUD6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mindin Polyclonal Antibody detects endogenous levels of Mindin protein.
Storage and Stability:	-20°C/1 year