

Antibodies Product Datasheet

Product Name:	JIK Antibody
Catalog Number:	C20696
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAO3. AA range:531-580
Concentration:	1 mg/ml
Gene Name:	TAOK3
Protein Name:	Serine/threonine-protein kinase TAO3
Other Name:	TAOK3; DPK; JIK; KDS; MAP3K18; Serine/threonine-protein kinase TAO3; Cutaneous T-cell lymphoma-associated antigen HD-CL-09; CTCL-associated antigen HD-CL-09; Dendritic cell-derived protein kinase; JNK/SAPK-inhibitory kinase; Jun kinase-inhi
Human Swiss Prot No.:	Q9H2K8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JIK Polyclonal Antibody detects endogenous levels of JIK protein.
Storage and Stability:	-20°C/1 year