

Antibodies Product Datasheet

Product Name:	DGK- β Antibody
Catalog Number:	C20741
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DGKB. AA range:657-706
Concentration:	1 mg/ml
Gene Name:	DGKB
Protein Name:	Diacylglycerol kinase beta
Other Name:	DGKB; DAGK2; KIAA0718; Diacylglycerol kinase beta; DAG kinase beta; 90 kDa diacylglycerol kinase; Diglyceride kinase beta; DGK-beta
Human Swiss Prot No.:	Q9Y6T7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DGK- β Polyclonal Antibody detects endogenous levels of DGK- β protein.
Storage and Stability:	-20°C/1 year