

Antibodies Product Datasheet

Product Name:	MKP-5 Antibody
Catalog Number:	C20742
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MKP-5 . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	DUSP10
Protein Name:	Dual specificity protein phosphatase 10
Other Name:	DUSP10; MKP5; Dual specificity protein phosphatase 10; Mitogen-activated protein kinase phosphatase 5; MAP kinase phosphatase 5; MKP-5
Human Swiss Prot No.:	Q9Y6W6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MKP-5 Polyclonal Antibody detects endogenous levels of MKP-5 protein.
Storage and Stability:	-20°C/1 year