

Antibodies Product Datasheet

Product Name:	FGF-BP2 Antibody
Catalog Number:	C20778
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FGFBP2. AA range:170-219
Concentration:	1 mg/ml
Gene Name:	FGFBP2
Protein Name:	Fibroblast growth factor-binding protein 2
Other Name:	FGFBP2; KSP37; Fibroblast growth factor-binding protein 2; FGF-BP2; FGF-binding protein 2; FGFBP-2; 37 kDa killer-specific secretory protein; Ksp37; HBp17-related protein; HBp17-RP
Human Swiss Prot No.:	Q9BYJ0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-BP2 Polyclonal Antibody detects endogenous levels of FGF-BP2 protein.
Storage and Stability:	-20°C/1 year