

Antibodies Product Datasheet

Product Name:	FGFR-5 Antibody
Catalog Number:	C20779
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from FGFR-5 . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	FGFRL1
Protein Name:	Fibroblast growth factor receptor-like 1
Other Name:	FGFRL1; FGFR5; FHFR; Fibroblast growth factor receptor-like 1; FGF receptor-like protein 1; FGF homologous factor receptor; FGFR-like protein; Fibroblast growth factor receptor 5; FGFR-5
Human Swiss Prot No.:	Q8N441
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGFR-5 Polyclonal Antibody detects endogenous levels of FGFR-5 protein.
Storage and Stability:	-20°C/1 year