

Antibodies Product Datasheet

Product Name:	HSP77/76 Antibody
Catalog Number:	C20789
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human HSP77/76.
Concentration:	1 mg/ml
Gene Name:	HSPA6/HSPA7
Protein Name:	Heat shock 70 kDa protein 6/Putative heat shock 70 kDa protein 7
Other Name:	HSPA7; HSP70B; Putative heat shock 70 kDa protein 7; Heat shock 70 kDa protein B; HSPA6; HSP70B'; Heat shock 70 kDa protein 6; Heat shock 70 kDa protein B'
Human Swiss Prot No.:	P48741
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSP77/76 Polyclonal Antibody detects endogenous levels of HSP77/76 protein.
Storage and Stability:	-20°C/1 year