

Antibodies Product Datasheet

Product Name:	IL-28R α Antibody
Catalog Number:	C20805
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human IL-28R α .
Concentration:	1 mg/ml
Gene Name:	IL28RA
Protein Name:	Interleukin-28 receptor subunit alpha
Other Name:	IL28RA; IFNLR1; LICR2; Interleukin-28 receptor subunit alpha; IL-28 receptor subunit alpha; IL-28R-alpha; IL-28RA; Cytokine receptor class-II member 12; Cytokine receptor family 2 member 12; CRF2-12; Interferon lambda receptor 1; IFN-lambda receptor 1; IFN
Human Swiss Prot No.:	Q8IU57
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-28R α Polyclonal Antibody detects endogenous levels of IL-28R α protein.
Storage and Stability:	-20°C/1 year