

Antibodies Product Datasheet

Product Name:	IFIT-5 Antibody
Catalog Number:	C20809
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IFIT5. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	IFIT5
Protein Name:	Interferon-induced protein with tetratricopeptide repeats 5
Other Name:	IFIT5; ISG58; RI58; Interferon-induced protein with tetratricopeptide repeats 5; IFIT-5; Interferon-induced 58 kDa protein; Retinoic acid- and interferon-inducible 58 kDa protein; P58
Human Swiss Prot No.:	Q13325
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFIT-5 Polyclonal Antibody detects endogenous levels of IFIT-5 protein.
Storage and Stability:	-20°C/1 year