

Antibodies Product Datasheet

Product Name:	IL-4 Antibody
Catalog Number:	C20839
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human IL-4.
Concentration:	1 mg/ml
Gene Name:	IL4
Protein Name:	Interleukin-4
Other Name:	IL4; Interleukin-4; IL-4; B-cell stimulatory factor 1; BSF-1; Binetrakin; Lymphocyte stimulatory factor 1; Pitakinra
Human Swiss Prot No.:	P05112
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-4 Polyclonal Antibody detects endogenous levels of IL-4 protein.
Storage and Stability:	-20°C/1 year