

Antibodies Product Datasheet

Product Name:	JAM-B Antibody
Catalog Number:	C20873
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human JAM-B.
Concentration:	1 mg/ml
Gene Name:	JAM2
Protein Name:	Junctional adhesion molecule B
Other Name:	JAM2; C21orf43; VEJAM; Junctional adhesion molecule B; JAM-B; Junctional adhesion molecule 2; JAM-2; Vascular endothelial junction-associated molecule; VE-JAM; CD322
Human Swiss Prot No.:	P57087
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	JAM-B Polyclonal Antibody detects endogenous levels of JAM-B protein.
Storage and Stability:	-20°C/1 year