

Antibodies Product Datasheet

Product Name:	Keratin 37/38 Antibody
Catalog Number:	C20876
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KRT37/38. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	KRT37/KRT38
Protein Name:	Keratin type I cuticular Ha7/Keratin type I cuticular Ha8
Other Name:	KRT37; HHA7; HKA7; KRTHA7; Keratin; type I cuticular Ha7; Hair keratin, type I Ha7; Keratin-37; K37; KRT38; HHA8; HKA8; KRTHA8; Keratin, type I cuticular Ha8; Hair keratin, type I Ha8; Keratin-38; K38
Human Swiss Prot No.:	O76014/O76015
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Keratin 37/38 Polyclonal Antibody detects endogenous levels of Keratin 37/38 protein.
Storage and Stability:	-20°C/1 year