

Antibodies Product Datasheet

Product Name:	Keratin 77 Antibody
Catalog Number:	C20879
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), WB
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Keratin 77
Concentration:	1 mg/ml
Gene Name:	KRT77 KRT1B
Protein Name:	Keratin 77
Other Name:	Keratin, type II cytoskeletal 1b (Cytokeratin-1B) (CK-1B) (Keratin-77) (K77) (Type-II keratin Kb39)
Human Swiss Prot No.:	Q7Z794
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Keratin 77
Storage and Stability:	-20°C/1 year