

Antibodies Product Datasheet

Product Name:	MYL3 Antibody
Catalog Number:	C20888
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYL3. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	MYL3
Protein Name:	Myosin light chain 3
Other Name:	MYL3; Myosin light chain 3; Cardiac myosin light chain 1; CMLC1; Myosin light chain 1; slow-twitch muscle B/ventricular isoform; MLC1SB; Ventricular/slow twitch myosin alkali light chain
Human Swiss Prot No.:	P08590
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MYL3 Polyclonal Antibody detects endogenous levels of MYL3 protein.
Storage and Stability:	-20°C/1 year