

Antibodies Product Datasheet

Product Name:	Group VI iPLA2 Antibody
Catalog Number:	C20892
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Group VI iPLA2.
Concentration:	1 mg/ml
Gene Name:	PLA2G6
Protein Name:	85/88 kDa calcium-independent phospholipase A2
Other Name:	PLA2G6; PLPLA9; 85/88 kDa calcium-independent phospholipase A2; Cal-PLA2; Group VI phospholipase A2; GVI PLA2; Intracellular membrane-associated calcium-independent phospholipase A2 beta; iPLA2-beta; Patatin-like phospholipase domain-contai
Human Swiss Prot No.:	O60733
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Group VI iPLA2 Polyclonal Antibody detects endogenous levels of Group VI iPLA2 protein.
Storage and Stability:	-20°C/1 year