

Antibodies Product Datasheet

Product Name:	KMT1B Antibody
Catalog Number:	C20905
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SUV39H2. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	SUV39H2
Protein Name:	Histone-lysine N-methyltransferase SUV39H2
Other Name:	SUV39H2; KMT1B; Histone-lysine N-methyltransferase SUV39H2; Histone H3-K9 methyltransferase 2; H3-K9-HMTase 2; Lysine N-methyltransferase 1B; Suppressor of variegation 3-9 homolog 2; Su(var)3-9 homolog 2
Human Swiss Prot No.:	Q9H5I1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KMT1B Polyclonal Antibody detects endogenous levels of KMT1B protein.
Storage and Stability:	-20°C/1 year