

Antibodies Product Datasheet

Product Name:	CSRNP2 Antibody
Catalog Number:	C20906
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAIP-12. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	CSRNP2
Protein Name:	Cysteine/serine-rich nuclear protein 2
Other Name:	CSRNP2; C12orf22; FAM130A1; TAIP12; Cysteine/serine-rich nuclear protein 2; CSRNP-2; Protein FAM130A1; TGF-beta-induced apoptosis protein 12; TAIP-12
Human Swiss Prot No.:	Q9H175
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CSRNP2 Polyclonal Antibody detects endogenous levels of CSRNP2 protein.
Storage and Stability:	-20°C/1 year