

Antibodies Product Datasheet

Product Name:	Cdc14a phosphatase Antibody
Catalog Number:	C20925
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Cdc14a phosphatase.
Concentration:	1 mg/ml
Gene Name:	CDC14A
Protein Name:	Dual specificity protein phosphatase CDC14A
Other Name:	CDC14A; Dual specificity protein phosphatase CDC14A; CDC14 cell division cycle 14 homolog A
Human Swiss Prot No.:	Q9UNH5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc14a phosphatase Polyclonal Antibody detects endogenous levels of Cdc14a phosphatase protein.
Storage and Stability:	-20°C/1 year