

Antibodies Product Datasheet

Product Name:	Chfr Antibody
Catalog Number:	C20927
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Chfr.
Concentration:	1 mg/ml
Gene Name:	CHFR
Protein Name:	E3 ubiquitin-protein ligase CHFR
Other Name:	CHFR; RNF196; E3 ubiquitin-protein ligase CHFR; Checkpoint with forkhead and RING finger domains protein; RING finger protein 196
Human Swiss Prot No.:	Q96EP1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Chfr Polyclonal Antibody detects endogenous levels of Chfr protein.
Storage and Stability:	-20°C/1 year