

Antibodies Product Datasheet

Product Name:	Cyclin M2 Antibody
Catalog Number:	C20929
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNNM2. AA range:571-620
Concentration:	1 mg/ml
Gene Name:	CNNM2
Protein Name:	Metal transporter CNNM2
Other Name:	CNNM2; ACDP2; Metal transporter CNNM2; Ancient conserved domain-containing protein 2; Cyclin-M2
Human Swiss Prot No.:	Q9H8M5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cyclin M2 Polyclonal Antibody detects endogenous levels of Cyclin M2 protein.
Storage and Stability:	-20°C/1 year