

Antibodies Product Datasheet

Product Name:	RFC2 Antibody
Catalog Number:	C20934
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RFC2. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	RFC2
Protein Name:	Replication factor C subunit 2
Other Name:	RFC2; Replication factor C subunit 2; Activator 1 40 kDa subunit; A1 40 kDa subunit; Activator 1 subunit 2; Replication factor C 40 kDa subunit; RF-C 40 kDa subunit; RFC40
Human Swiss Prot No.:	P35250
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RFC2 Polyclonal Antibody detects endogenous levels of RFC2 protein.
Storage and Stability:	-20°C/1 year