

Antibodies Product Datasheet

Product Name:	Bcl-7a Antibody
Catalog Number:	C20948
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BCL7A. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	BCL7A
Protein Name:	B-cell CLL/lymphoma 7 protein family member A
Other Name:	BCL7A; B-cell CLL/lymphoma 7 protein family member A
Human Swiss Prot No.:	Q4VC05
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bcl-7a Polyclonal Antibody detects endogenous levels of Bcl-7a protein.
Storage and Stability:	-20°C/1 year