

Antibodies Product Datasheet

Product Name:	KLF12 Antibody
Catalog Number:	C20956
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human KLF12.
Concentration:	1 mg/ml
Gene Name:	KLF12
Protein Name:	Krueppel-like factor 12
Other Name:	KLF12; AP2REP; HSPC122; Krueppel-like factor 12; Transcriptional repressor AP-2rep
Human Swiss Prot No.:	Q9Y4X4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KLF12 Polyclonal Antibody detects endogenous levels of KLF12 protein.
Storage and Stability:	-20°C/1 year