

Antibodies Product Datasheet

Product Name:	CYP4V2 Antibody
Catalog Number:	C20985
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CYP4V2.
Concentration:	1 mg/ml
Gene Name:	CYP4V2
Protein Name:	Cytochrome P450 4V2
Other Name:	CYP4V2; Cytochrome P450 4V2
Human Swiss Prot No.:	Q6ZWL3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP4V2 Polyclonal Antibody detects endogenous levels of CYP4V2 protein.
Storage and Stability:	-20°C/1 year