

Antibodies Product Datasheet

Product Name:	CYP4Z1/2 Antibody
Catalog Number:	C21021
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CYP4Z1/2.
Concentration:	1 mg/ml
Gene Name:	CYP4Z1/CYP4Z2P
Protein Name:	Cytochrome P450 4Z1/Putative cytochrome P450 family member 4Z2
Other Name:	
Human Swiss Prot No.:	Q5VVE4/Q8N1L4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP4Z1/2 Polyclonal Antibody detects endogenous levels of CYP4Z1/2 protein.
Storage and Stability:	-20°C/1 year