

Antibodies Product Datasheet

Product Name:	CYP2A6/7/13 Antibody
Catalog Number:	C21054
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CYP2A6/7/13.
Concentration:	1 mg/ml
Gene Name:	CYP2A6/CYP2A7/CYP2A13
Protein Name:	Cytochrome P450 2A6/Cytochrome P450 2A7/Cytochrome P450 2A13
Other Name:	CYP2A6; CYP2A3; Cytochrome P450 2A6; 1,4-cineole 2-exo-monooxygenase; CYP2A6; Coumarin 7-hydroxylase; Cytochrome P450 IIA3; Cytochrome P450(I); CYP2A7; Cytochrome P450 2A7; CYP2A7; Cytochrome P450 IIA4; CYP2A13; Cytochrome P450 2A13; CYP2A13
Human Swiss Prot No.:	P11509/P20853/Q16696
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP2A6/7/13 Polyclonal Antibody detects endogenous levels of CYP2A6/7/13 protein.
Storage and Stability:	-20°C/1 year