

Antibodies Product Datasheet

Product Name:	SIRT2 Antibody
Catalog Number:	C21067
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SIRT2. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	SIRT2
Protein Name:	NAD-dependent protein deacetylase sirtuin-2
Other Name:	SIRT2; SIR2L; SIR2L2; NAD-dependent protein deacetylase sirtuin-2; Regulatory protein SIR2 homolog 2; SIR2-like protein 2
Human Swiss Prot No.:	Q8IXJ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SIRT2 Polyclonal Antibody detects endogenous levels of SIRT2 protein.
Storage and Stability:	-20°C/1 year