

Antibodies Product Datasheet

Product Name:	KIR2.3 Antibody
Catalog Number:	C21080
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KCNJ4. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	KCNJ4
Protein Name:	Inward rectifier potassium channel 4
Other Name:	KCNJ4; IRK3; Inward rectifier potassium channel 4; HIRK2; HRK1; Hippocampal inward rectifier; HIR; Inward rectifier K(+) channel Kir2.3; IRK-3; Potassium channel; inwardly rectifying subfamily J member 4
Human Swiss Prot No.:	P48050
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KIR2.3 Polyclonal Antibody detects endogenous levels of KIR2.3 protein.
Storage and Stability:	-20°C/1 year