

Antibodies Product Datasheet

Product Name:	Bag-4 Antibody
Catalog Number:	C21104
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human Bag-4.
Concentration:	1 mg/ml
Gene Name:	BAG4
Protein Name:	BAG family molecular chaperone regulator 4
Other Name:	BAG4; SODD; BAG family molecular chaperone regulator 4; BAG-4; Bcl-2-associated athanogene 4; Silencer of death domains
Human Swiss Prot No.:	O95429
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bag-4 Polyclonal Antibody detects endogenous levels of Bag-4 protein.
Storage and Stability:	-20°C/1 year