

Antibodies Product Datasheet

Product Name:	Brn-3a Antibody
Catalog Number:	C21106
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human POU4F1. AA range:311-360
Concentration:	1 mg/ml
Gene Name:	POU4F1
Protein Name:	POU domain class 4 transcription factor 1
Other Name:	POU4F1; BRN3A; RDC1; POU domain; class 4, transcription factor 1; Brain-specific homeobox/POU domain protein 3A; Brain-3A; Brn-3A; Homeobox/POU domain protein RDC-1; Oct-T1
Human Swiss Prot No.:	Q01851
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brn-3a Polyclonal Antibody detects endogenous levels of Brn-3a protein.
Storage and Stability:	-20°C/1 year