

Antibodies Product Datasheet

Product Name:	MAGE-C2 Antibody
Catalog Number:	C21110
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from MAGE-C2 . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	MAGEC2
Protein Name:	Melanoma-associated antigen C2
Other Name:	MAGEC2; HCA587; MAGEE1; Melanoma-associated antigen C2; Cancer/testis antigen 10; CT10; Hepatocellular carcinoma-associated antigen 587; MAGE-C2 antigen; MAGE-E1 antigen
Human Swiss Prot No.:	Q9UBF1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAGE-C2 Polyclonal Antibody detects endogenous levels of MAGE-C2 protein.
Storage and Stability:	-20°C/1 year