

Antibodies Product Datasheet

Product Name:	Lsk Antibody
Catalog Number:	C21131
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Lsk.
Concentration:	1 mg/ml
Gene Name:	MATK
Protein Name:	Megakaryocyte-associated tyrosine-protein kinase
Other Name:	MATK; CTK; HYL; Megakaryocyte-associated tyrosine-protein kinase; CSK homologous kinase; CHK; Hematopoietic consensus tyrosine-lacking kinase; Protein kinase HYL; Tyrosine-protein kinase CTK
Human Swiss Prot No.:	P42679
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Lsk Polyclonal Antibody detects endogenous levels of Lsk protein.
Storage and Stability:	-20°C/1 year