

Antibodies Product Datasheet

Product Name:	DRAK2 Antibody
Catalog Number:	C21136
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human STK17B. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	STK17B
Protein Name:	Serine/threonine-protein kinase 17B
Other Name:	STK17B; DRAK2; Serine/threonine-protein kinase 17B; DAP kinase-related apoptosis-inducing protein kinase 2
Human Swiss Prot No.:	O94768
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DRAK2 Polyclonal Antibody detects endogenous levels of DRAK2 protein.
Storage and Stability:	-20°C/1 year