

Antibodies Product Datasheet

Product Name:	HUNK Antibody
Catalog Number:	C21146
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HUNK. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	HUNK
Protein Name:	Hormonally up-regulated neu tumor-associated kinase
Other Name:	HUNK; MAKV; Hormonally up-regulated neu tumor-associated kinase; B19; Serine/threonine-protein kinase MAK-V
Human Swiss Prot No.:	P57058
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HUNK Polyclonal Antibody detects endogenous levels of HUNK protein.
Storage and Stability:	-20°C/1 year