

Antibodies Product Datasheet

Product Name:	MASTL Antibody
Catalog Number:	C21155
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human MASTL.
Concentration:	1 mg/ml
Gene Name:	MASTL
Protein Name:	Serine/threonine-protein kinase greatwall
Other Name:	MASTL; GW; GWL; THC2; Serine/threonine-protein kinase greatwall; GW; GWL; hGWL; Microtubule-associated serine/threonine-protein kinase-like; MAST-L
Human Swiss Prot No.:	Q96GX5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MASTL Polyclonal Antibody detects endogenous levels of MASTL protein.
Storage and Stability:	-20°C/1 year