

Antibodies Product Datasheet

Product Name:	MerTK Antibody
Catalog Number:	C21157
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MERTK. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	MERTK
Protein Name:	Tyrosine-protein kinase Mer
Other Name:	MERTK; MER; Tyrosine-protein kinase Mer; Proto-oncogene c-Mer; Receptor tyrosine kinase MerTK
Human Swiss Prot No.:	Q12866
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MerTK Polyclonal Antibody detects endogenous levels of MerTK protein.
Storage and Stability:	-20°C/1 year