

Antibodies Product Datasheet

Product Name:	MRCK α Antibody
Catalog Number:	C21158
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human MRCK α .
Concentration:	1 mg/ml
Gene Name:	CDC42BPA
Protein Name:	Serine/threonine-protein kinase MRCK alpha
Other Name:	CDC42BPA; KIAA0451; Serine/threonine-protein kinase MRCK alpha; CDC42-binding protein kinase alpha; DMPK-like alpha; Myotonic dystrophy kinase-related CDC42-binding kinase alpha; MRCK alpha; Myotonic dystrophy protein kinase-like alpha
Human Swiss Prot No.:	Q5VT25
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRCK α Polyclonal Antibody detects endogenous levels of MRCK α protein.
Storage and Stability:	-20°C/1 year