

Antibodies Product Datasheet

Product Name:	Duo Antibody
Catalog Number:	C21201
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Duo . at AA range: 810-890
Concentration:	1 mg/ml
Gene Name:	KALRN
Protein Name:	Kalirin
Other Name:	KALRN; DUET; DUO; HAPIP; TRAD; Kalirin; Huntingtin-associated protein-interacting protein; Protein Duo; Serine/threonine-protein kinase with Dbl- and pleckstrin homology domain
Human Swiss Prot No.:	O60229
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Duo Polyclonal Antibody detects endogenous levels of Duo protein.
Storage and Stability:	-20°C/1 year