

Antibodies Product Datasheet

Product Name:	CaMKK2 Antibody
Catalog Number:	C21217
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CAMKK2. AA range:381-430
Concentration:	1 mg/ml
Gene Name:	CAMKK2
Protein Name:	Calcium/calmodulin-dependent protein kinase kinase 2
Other Name:	CAMKK2; CAMKKB; KIAA0787; Calcium/calmodulin-dependent protein kinase kinase 2; CaM-KK 2; CaM-kinase kinase 2; CaMKK 2; Calcium/calmodulin-dependent protein kinase kinase beta; CaM-KK beta; CaM-kinase kinase beta; CaMKK beta
Human Swiss Prot No.:	Q96RR4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CaMKII β Polyclonal Antibody detects endogenous levels of CaMKII β protein.
Storage and Stability:	-20°C/1 year